



(407) 614-8781 office
info@PaladinHC.com

CHALLENGE:

To deliver Oxygen to multiple patients in unconventional environments due to catastrophe or mass casualty events.



www.paladinhc.com

MASS CASUALTY MEDICAL GAS DELIVERY SYSTEM



The safe and effective delivery of Oxygen for mass casualty or major fire incidents has always been a challenge.

Paladin's Surge Gas Delivery System meets those needs. It deploys in minutes, sets up in rough terrain, and can simultaneously supply 8 victims with Oxygen. It stores in a portable and easy to carry case.

Kit Consists of:

- ☒ (2) Tripods with variable height legs for uneven terrain
- ☒ (1) Collapsible 8' Paladin Evolution Medical Equipment Rail (breaks down to 4')
- ☒ (7) Rail mounted oxygen blocks with flow meters and DISS check valve
- ☒ (7) Hoses to connect each block and to main O₂ source
- ☒ (1) Carrying case for easy storage and deployment

Paladin Healthcare LLC, Oakland, FL (407) 614-8781